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Oblique View of Petri Dish

Side View of Petri Dish



Figure 2. Some culture plates needed for IVP. From left to right are a Falcon 60 x 15 petri dish used for preparing microdrops, a Nunclon 4 well plate for fertilization, a Fisher X-plate used in several steps for cleaning up preparations of COCs and embryos, and a Fisher intergrid plate for searching for cumulus-oocyte complexes. All can be purchased from [Fisher](#).

**Figure 3.** Preparation of microdrops. Shown here are 9-50  $\mu$ l drops on the bottom of a 60 mm petri dish. To complete the drops, they are covered with mineral oil by depositing oil using a Pasteur pipette.


**Figure 4.** Some items needed for oocyte collection. Shown on the warming plate are 2 1-L plastic beakers (one containing cleaned ovaries and one to collect ovaries after processing) and a 400 ml beaker containing ~100 ml OCM. Also shown are a hemostat, scalpel handle, disposable scalpel blades and bench paper.

**Figure 5.** *Massaging the ovary to remove blood.*


**Figure 6.** Slashing of ovaries to recover oocytes. In panel A, a hemostat is attached to the base of the ovary to hold the ovary firmly in place. The excess tissue from the ovarian stalk is removed in panels


© Rivera et al., 2000

**Figure 7.** Harvesting of oocytes from slashed ovaries. In the left panel, the slashed ovary is being swirled in OCM. In the right panel, the ovary is being pressed against the side of the beaker to allow drainage of OCM.


**Figure 8.** Cleaning up the preparation of oocytes collected by ovary slashing. Panel A depicts the placement of tubes in the water bath for settling and Panel B the supplies necessary for this step. Aspiration of oocytes from the bottom of the tube is shown in Panel C while Panel D represents the technique used to pour the aspirated pellet through the cell strainer. Rinsing of oocytes from the cell strainer into an integrid petri dish is shown in Panel E.


**Figure 9.** Searching for oocytes. After COCs have been washed, the contents are poured into an Intergrid plate and examined under a dissecting scope for COCs. These can be picked up using a wiretrol or other device.

Figure 10. An X-plate used to wash COCs, purchased from Fisher.

**Figure 12.** Some instruments used to pick up embryos and oocytes. From left to right are 1) a 1 cc syringe with an extension of rubber tubing connected to a [Unopette](#), 2) a wiretrol (from Drummond Scientific; we purchase from [Fisher](#)), 3) the same device as #1 without the rubber tubing extension and 4) a 5 ml Drummond Microdispenser (also purchased from [Fisher](#)). For tips on how to use these instruments, [click here](#).

**Figure 13.** Transfer of COCs to microdrops.

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




**Fig. 2.11.** Effect of recovery method on the diameters of good-quality oocytes obtained (from Carolan et al., 1994).

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**FIGURE 6** Diagram of transvaginal ultrasound-guided follicular aspiration in cattle. Ovary (1) is positioned against the transducer by rectal manipulation (2) while the ultrasound transducer containing the needle guide and aspiration needle (4) is placed against the anterior portion of the vagina beside the cervix (3). (Kruip *et al.*, 1991. *Vet. Rec.* **128**, 208–210.)

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**Table 5.3.** Composition of media M-199, KRB, mKRB and CZB

Components	Medium <sup>a</sup>			
	M-199 <sup>b</sup>	KRB	mKRB	CZB
NaCl	8.000	6.721	6.721	4.769
KCl	0.400	0.356	0.356	0.360
NaHCO <sub>3</sub>	0.350	2.100	2.106	2.110
KH <sub>2</sub> PO <sub>4</sub>	0.060	0.162	0.162	0.161
MgSO <sub>4</sub> ·7H <sub>2</sub> O	—	0.293	0.294	0.294
Na <sub>2</sub> HPO <sub>4</sub>	0.048	—	—	—
CaCl <sub>2</sub>	0.140	0.347	0.189	—
CaCl <sub>2</sub> ·2H <sub>2</sub> O	—	—	—	0.249
Na pyruvate	—	—	—	0.029
BSA	varies	—	4.000	5.000
Na lactate	—	—	—	3.508
Glucose	1.000	—	1.000	—
EDTA	—	—	—	0.041
Glutamine	—	—	—	1.000

<sup>a</sup>Components are in g l<sup>-1</sup> except for glutamine, which is added at a concentration of 1.0 mM. Penicillin and streptomycin were added at concentrations of 100 iu ml<sup>-1</sup> and 50 μg ml<sup>-1</sup>, respectively.

<sup>b</sup>M-199 also contains certain amino acids, vitamins and trace elements.

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**Table 4.8.** SOF formulation employed as a single culture system for bovine embryo production (from Gandhi *et al.*, 2000).

Reagent (mmol l <sup>-1</sup> )	SOFM (maturation)	SOFF (fertilization)	SOFS (sperm wash)	SOFC1 (culture 1)	SOFC2 (culture 2)
Sodium pyruvate	0.33	0.33	1.0	0.33	0.33
Glutamine	1.0	–	–	1.0	1.0
Glucose	1.5	–	–	1.5	3.0
Sodium lactate	–	–	18.3	–	–
HEPES	–	–	12.5	–	–
Bovine serum albumin	–	6 mg/ml fatty acid free	3 mg/ml fraction V	8 mg/ml crystallized	8 mg/ml crystallized
EDTA	–	–	–	0.1	–
Taurine	–	–	–	0.1	–
MEM. NEAA	1x	2x	–	1x	1x
MEM. EAA	1x	–	–	1x	1x
MEM. vitamins	–	–	–	–	1x

MEM, Eagle's minimum essential medium; NEAA, non-essential amino acids; EAA, essential amino acids. Supplements added to base SOF medium for specific developmental stages.

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
IVF JAPAN

## Indication and Procedure

Indications for IVM	PCO, Normal menstrual cycle, poor quality embryo case
Monitoring	From Day 7
Timing of oocyte collection	Two follicles >7mm
Cancellation criteria	Dominant follicle, Cyst formation
HCG priming	10000iu HCG 36 hrs before OPU
Oocyte collection	Aspiration pressure 100mmHG IVF OSAKA IVM needle
IVM medium	IVM Medium (Medicult), 10% Patient's serum, HCG 100 IU/ l, FSH 75 IU/l
Maturation time	24 hours
Fresh or Frozen Embryo transfer	Endom.>8mm Fresh ET

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
**Figure 3** Ultrasound image taken during oocyte retrieval. The highly echogenic band around the distal end of the needle and the tip of the follicle aspiration needle are visualized in the superior-most follicle. The needle is maneuvered within the ovary to aspirate all follicles. *Source:* Image courtesy of Dr. Roger Stronell.

**Figure 4** Midsagittal view of the uterus. The cervix is to the right of the image and the fundus is to the left. The endometrium is well demarcated and shows a pronounced, thick “triple-line” pattern associated with a higher probability of implantation following embryo transfer.

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














**Fig 9.11** After follicle fluid aspiration, a cumulus–corona–oocyte complex (COC) has to be found floating between other granulosa cell debris (arrow). Arrow in (a), (b), and (c) shows the COC from low to large magnification as visualized under the stereomicroscope. Bar is 500 µm.

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

**Fig 9.12** Classification of cumulus–corona–oocyte complex (COC) retrieved from small antral follicles. (a) CM1, CEO, CO2: the oocyte is surrounded by a compact mass of four to five layers of granulosa cells. The germinal vesicle is clearly visible within the oocyte at a peripheral position. (b) CM1, CE2, CO2: there is expansion of the distal layer of granulosa cells (cumulus), the proximal granulosa cells surrounding the oocyte are still compacted. (c) CM1, CE1, CO2: there is expansion of both distal (cumulus) and proximal (corona cells) layers of granulosa cells. (d) Distal granulosa cells are dispersed. The oocyte is partially naked and connections between cumulus cells and oocytes are lost.

**Fig 9.13** Photomicrograph of a cumulus–oocyte complex with an expanded cell morphology pattern (CM2, CE2, CO2). (a) Nuclear maturity of oocyte is not visible under an inverted microscope. (b) After cell ‘spreading,’ the GV is visualized within the ooplasm (arrow).

**Fig 9.14** Spindle images of human in vitro matured oocytes using the PolScope software showing different locations of spindle in the ooplasm in relation to polar body position. (Courtesy of Lin Liu and David Keefe, University of Florida, Florida.)

**Table 7.1 Oocyte–corona–cumulus complex evaluation scheme**

Groups	OCCC morphology
Mature	Expanded cumulus Radiant corona Distinct zona pellucida, clear ooplasm Expanded well-aggregated membrana granulosa cells
Approximately mature	Expanded cumulus mass Slightly compact corona radiata Expanded well-aggregated membrana granulosa cells
Immature	Dense compact cumulus if present Adherent compact layer of corona cells Ooplasm if visible with the presence of germinal vesicle Compact and non-aggregated membrana granulosa cells
Post-mature	Expanded cumulus with clumps Radiant corona radiata, yet often clumped, irregular, or incomplete Visible zona, slightly granular or dark ooplasm Small and relatively nonaggregated membrana granulosa cells
Atretic	Rarely with associated cumulus mass Clumped and very irregular corona radiata if present Visible zona, dark and frequently misshapen ooplasm Membrana granulosa cells with very small clumps of cells

Adapted from Lin et al.<sup>95</sup>

**Table 9.5 Maturity grading of cumulus–oocyte complexes (COCs) and oocyte (after cumulus–corona denudation) during prematuration days in culture**

1. *Grading of granulosa cell mass: cumulus expansion and oocyte coverage*

Cumulus mass (CM):

3 or fewer layers of cumulus cells	(CM0)
more than 3 but fewer than 10 layers of cumulus cells	(CM1)
10 or more layers of cumulus cells	(CM2)

Cumulus expansion (CE):

tight, dense cells	(CE0)
moderate expansion of cells	(CE1)
fully expanded cells	(CE2)

Contact (CO) between cumulus cells and oocyte:

naked	(CO0)
partially naked	(CO1)
fully enclosed	(CO2)
fully enclosed and part of follicle wall	(CO3)

2. *Assessment of oocyte nuclear maturation stage*

GV

GVBD

PB

3. *Assessment of oocyte morphology*

Oocyte diameter (µm)

Oocyte cytoplasm

presence of inclusions:

vacuoles/refractile bodies

darkness:

clear/dark

granularity:

homogeneous/granular

Zona:

normal/abnormal

Perivitelline space:

normal/enlarged

Oocyte shape:

regular/irregular

Polar body:

intact/fragmented

GV, germinal vesicle; GVBD, germinal vesicle breakdown; PB, polar body

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**Fig 7.8** Giant oocyte (arrow) and normally sized oocyte. (Magnification 400x.)

**Fig 7.9** Metaphase II oocytes obtained after ovarian hyperstimulation, all with approximately the same morphological appearance. (Magnification 400x.)

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Vial 1: The LAG Medium is for preincubation of immature oocytes

Vial 2: The IVM® Medium is a basal medium for maturing immature oocytes with supplements

MediCult IVM® is for use

- in healthy regularly cycling women referred due to male factor and/or tubal disease
- women with PCO

## Composition

Vial 1, LAG Medium:

- Synthetic Serum Replacement (SSR®) (USA: ART Supplement)
- Physiological salts
- Human serum albumin (HSA)
- Glucose
- Sodium pyruvate
- Sodium bicarbonate
- Penicillin
- Streptomycin
- Phenol Red

Vial 2, IVM® Medium:

- Glucose
- Sodium Pyruvate
- Sodium Acetate
- Cholesterol
- Physiological salts
- Amino acids
- Nucleotides
- Vitamins
- Sodium bicarbonate
- Penicillin
- Streptomycin
- Phenol Red

IVF JAPAN

## Medium for IVM -Medicult

- Vial 1(Lag medium)  
for preincubation
- Vial 2(IVM Medium)  
a basal medium for maturing immature  
oocyte with supplements


IVF JAPAN

## Medium for IVM -Medicult

- Vial 2, IVM Medium
  - Glucose and derived metabolism
  - Physiological salts
  - Essential amino acids
  - Non-essential amino acid
  - Vitamins
  - DNA bases
  - Sodium pyruvate
  - Sodium bicarbonate
  - Streptomycin 50mg/l
  - Penicillin 50000IU/l
  - Phenol Red

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
ISM  
42-44 h


## TYPES OF IVM

**Fig 9.10** Diagram representing the diversified IVM clinical protocols applied in PCO, PCOS, or normo-ovulatory women. These protocols reflect on morphology and maturation timing of oocytes retrieved. (a) Immature oocytes are retrieved prior to the formation of a dominant follicle during a natural cycle: oocytes retrieved are enclosed in a compact mass of cumulus cells and the majority of the oocytes reach maturation between 28 h and 36 hours of culture. (b) Immature oocytes are collected following a mild stimulation: from patients primed only with FSH, the oocytes retrieved have a full compacted cell morphology and the majority of oocytes mature within 28–36 hours of culture; with the administration of hCG with or without FSH priming, maturation is more rapid, within 24–30 hours of culture, and the oocytes retrieved from the most advanced follicles have a more expanded cell morphology. (c) Immature oocytes collected from patients undergoing controlled ovarian hyperstimulation for IVF/ICSI treatment are often in vitro matured in the absence of adjacent surrounding granulosa cells and the majority of the oocytes mature within a time frame of 24–30 hours.

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IVF JAPAN

# Protocol of IVM

Oocyte collection **ICSI**

Embryo transfer

↑  
Ultrasonogram

Estrogen patch 6 mg every other day

HCG 10,000 I.U.(36 hrs prior to OPU)

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