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dr n wet R. Faúndez

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
Sperm DNA fragmentation: mechanisms of origin, impact on reproductive outcome, and analysis

Denny Sakkas, Ph.D., and Juan G. Alvarez, M.D., Ph.D, Fertility and Sterility Vol. 93, No. 4, March 1, 2010

FIGURE 1

Major mechanisms of inducing DNA damage in spermatozoa during either the production or the transport of sperm cells: (i) apoptosis during the process of spermatogenesis; (ii) DNA strand breaks produced during the remodelling of sperm chromatin during the process of spermiogenesis; (iii) post-testicular DNA fragmentation induced, mainly by oxygen radicals, during sperm transport through the seminiferous tubules and the epididymis (increasing DNA damage is indicated by size of red flashes and gradient darkening in tract); (iv) DNA fragmentation induced by endogenous caspases and endonucleases; (v) DNA damage induced by radiotherapy and chemotherapy; and (vi) DNA damage induced by environmental toxicants.

Sakkas. Sperm DNA fragmentation. Fertil Steril 2010.

TABLE 1

Outcome results of TESA-ICSI in couples with repeated IVF failure and high levels of sperm DNA fragmentation in semen.

Characteristic	Value
Semen and case parameters	
Cases evaluated	68
Cases of TESA-ICSI performed	31
Mean DFI in semen (%)	39.4
Mean sperm concentration (million/mL)	44.7
Mean male age (y)	41.9
Outcome results	
Fertilization rate (%)	58 (20-100)
Number of embryos transferred	2.3
Clinical pregnancy rate (%)	40.0
Pregnancy rate in the first TESA-ICSI cycle (%)	93

Note: DFI = DNA fragmentation index; TESA-ICSI = testicular sperm extraction-intracytoplasmic sperm injection.

Sakkas. Sperm DNA fragmentation. Fertil Steril 2010.

TABLE 2

Pregnancy outcome in patients with high levels of sperm DNA fragmentation when treated with ejaculated compared with testicular sperm.

Sperm	Biochemical pregnancy rate (%)	Clinical pregnancy rate (%)	Implantation rate (%)	Miscarriage rate (%)
Ejaculated (n = 42)	6.90	13.79	6.56	75.0
Testicular (n = 30)	2.5	40.0	28.09	6.25
P value	NS	.035	.0021	.017

Note: All patients included in the study had a DNA fragmentation index value in semen by TUNEL >20%; statistical analysis by χ^2 test.

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TEST TUNEL

(TERMINAL DEOKSYNUCLEOTIDYL TRANSFERASES dUTP AND LABELING)

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



TEST COMET

(SCGE - SINGLE CELL GEL ELECTROPHORESIS)

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Alkaline Lysis


Alkaline Unwinding


Single Strand
and
Double Strand
Break Detection

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

Obraz DNA:
50% DNA w ogonie











DNA BREAK DETECTION FLUORESCENT *IN SITU* HYBRIDIZATION TEST (DBD FISH)

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application


keep dry

Halosperm® Kit from Halotech Dna, SL is a simple test that allows assessment of sperm DNA fragmentation in humans.

Principle of the method:

The method is based on the Sperm Chromatin Dispersion (SCD) test (Fernández et al., J. Androl 24: 59-66, 2003; Fertil Steril 84: 833-842, 2005). Intact unfixed spermatozoa (fresh, frozen/unthawed, diluted samples) are immersed in an inert agarose microgel on a pretreated slide. An initial acid treatment denatures DNA in those sperm cells with fragmented DNA. Following this, the lysing solution removes most of the nuclear proteins, and in the absence of massive DNA breakage produces nucleoids with large halos of spreading DNA loops, emerging from a central core. However, the nucleoids from spermatozoa with fragmented DNA either do not show a dispersion halo or the halo is minimal.


Fig. 10.3 High magnification selected sperm (a, b) and the characterization of SDF in the same sperm (a', b'). The SCD test allows the direct assessment of the DNA status and the sperm morphology

Comparison of sperm DNA fragmentation tests



In order to evaluate sperm DNA fragmentation, Halotech has developed the Halosperm[®] kit: a fast and simple in vitro diagnostic kit based on the Sperm Chromatin Dispersion (SCD) technique first described by Fernández et al., 2003. Here is how our kit compares to other techniques available on the market:

Underlying principles of the main sperm DNA fragmentation tests

ASSAY	PRINCIPLE	MEASURED PARAMETER	DETECTION METHOD
TUNEL	Addition of labeled dUDP nucleotides with terminal deoxynucleotidyl transferase to both SS and DS DNA breaks. Template-independent	% cells with labeled DNA	Fluorescence or bright field microscopy / Flow cytometry
COMET	Electrophoresis of sperm cells reveals fragmented DNA. Alkaline conditions denature DNA to reveal SS and DS DNA breaks. Neutral conditions reveal mostly DS breaks	% cells with migration tails containing fragmented DNA (also length of tail, % of DNA in tail)	Fluorescence microscopy
SCSA	Mild acid treatment denatures DNA with DS or SS breaks. Acridine orange binds DNA and fluoresces green with DS DNA (non-denatured) and red with SS DNA (denatured)	DNA Fragmentation Index (DFI): cells with red fluorescence divided by total number of cells (red +green). Expressed as %	Flow cytometry
SCD	Mild acid denaturation of DNA and lysis of protamines creates a chromatin decondensation halo around sperm head when DNA is intact and no halo when DNA is damaged	% cells with small or absent decondensation halo	Bright field microscopy


TUNEL, TdT (Terminal deoxynucleotidyl Transferase) –mediated dUDP nick-end labeling. SCSA, Sperm Chromatin Structure Assay. SCD, Sperm Chromatin Decondensation test, SS, Single-Stranded. DS, Double-Stranded. (Tarozzi et al., 2007, Zini and Sigman, 2009)–

The results obtained with Halosperm® are comparable to those of TUNEL and SCSA

Sperm DNA fragmentation was analyzed in fertile donors ($n=7$, solid line) and infertile patients ($n=60$, dotted line). The correlations between the techniques (left, SCD and TUNEL, and right, SCD and SCSA) were analyzed using Student's t-test and are shown above each graph. (Adapted from: Chohan et al., 2006)

The same sperm cells are stained using TUNEL and Halosperm®

Sperm cells were embedded onto a slide and processed using the TUNEL assay followed by the SCD test. (A,B,C,D – left picture) Before Diff-quick staining, TUNEL positive sperms were stained brown. (A,B,C,D – right picture) After Diff-quick staining, TUNEL-positive sperms showed no halo and TUNEL-negative sperms showed a large halo. Pinhead sperm cells with no halo (Arrow) could not be observed with the TUNEL assay before Diff-quick staining. (Adapted from: Zhang et al., 2009).

Comparison of advantages and disadvantages for each technique

ASSAY	MAIN ADVANTAGES	MAIN DISADVANTAGES
TUNEL	<ol style="list-style-type: none"> 1. Can perform on few sperm 2. For certain protocols use of fluorescence microscopy or flow cytometry may not be necessary 	<ol style="list-style-type: none"> 1. Thresholds not standardized 2. Variable assay protocols 3. Not specifically designed for human spermatozoa 4. Labor intensive
COMET	<ol style="list-style-type: none"> 1. Sensitive 2. Can perform on few sperm 	<ol style="list-style-type: none"> 1. Requires complicated imaging software 2. Variable assay protocols 3. Unclear thresholds
SCSA	<ol style="list-style-type: none"> 1. Many cells rapidly examined 2. Most published data 	<ol style="list-style-type: none"> 1. Expensive equipment necessary 2. Small variations in lab conditions affect results 3. Calculations involve qualitative decisions

(Adapted from Zini and Sigman, 2009)

List of competitive advantages offered by SCD technique - Halosperm® kit

1. Has been specifically designed for use with human spermatozoa and unlike some other techniques, there exists only one standard protocol that is highly reproducible. For instance, it was recently demonstrated that small variations in crucial steps of the TUNEL assay can greatly affect measures of sperm DNA fragmentation (Muratori et al., 2009).
2. Is one of the most economical techniques to measure sperm DNA fragmentation on the market.
3. Is the fastest technique available on the market. The turnaround time is under 1 hour, and multiple slides can be processed at a time.
4. Is the simplest and easiest technique available on the market. Only basic laboratory equipment is needed to perform the test, and slides can be analyzed with any bright field light microscope.
5. Can be used with few sperm cells, such as in cases of oligozoospermia.
6. Is a highly versatile tool that can also be used in research in conjunction with a number of other techniques, e.g. labeling for 8-oxoguanine to measure oxidative stress, FISH to measure chromosomal aberrations or immunocytochemistry to investigate protein status,...

Kit for DNA fragmentation analysis in *Canis familiaris* spermatozoa (10 determinations)

KIT CONTAINS

- 10 precoated slides
- 10 vials with low melting agarose
- 1 flask with 100 ml of lysing solution. **Work under continuous air removal or cover trays during incubation step.**

BACKGROUND

The **SPERM-HALOMAX[®]** kit from **HALOTECH DNA SL** bursts into the market as a pioneering product for the analysis of fragmentation levels in spermatozoa both in domestic and experimental animals.

SPERM-HALOMAX[®] is based on the differential response of fragmented and unfragmented spermatozoa nuclei to a protein depletion treatment. The extraction of nuclear proteins from spermatozoa containing fragmented DNA releases DNA fragments between two breakage points. Sperm nuclei disperse chromatin around, forming a low stained peripheral halo, which is distinguishable under low magnification microscopy. These images contrast sharply with the intensely stained nuclei of the spermatozoa which occupies a centered position. On the contrary, those spermatozoa containing unfragmented DNA develop a very small dispersion halo, and appears as a very thin crown around the core. **SPERM HALOMAX[®]** allows the quick assessment of the proportion of spermatozoa with fragmented DNA.

MATERIAL AND EQUIPMENT REQUIRED BUT NOT PROVIDED

Fluorescence microscope, refrigerator at 4°C, water bath at 37°C and 90-100°C, latex gloves, glass coverslips (18x18 mm, 22x22mm, 24x60 mm), micropipettes, plastic tank for horizontal incubations, distilled water, 70%, 90% and 100% ethanol, staining solutions.

RECOMMENDED STAINING SOLUTIONS

For fluorescence microscopy: any DNA-specific binding fluorochrome could be used (DAPI, propidium iodide, ethidium bromide or synergy brand related fluorochromes). Antifading solution.

INSTRUCTIONS FOR USE

Including the sperm sample in agarose microgel

Set the lysing buffer at room temperature (22°C)

- 1) Dilute the sperm sample in culture medium or Phosphate Buffered Saline (PBS), to give a final concentration of 5 - 10 million spermatozoa per millilitre. Either fresh or liquid nitrogen frozen samples can be used.
- 2) In order to melt the agarose, place vial with low melting agarose onto a float and both into a water bath at 90-100°C for 5 minutes.

- 3) Transfer vial together with the float to a thermostatic water bath at 37°C, and leave it for 5 min to equilibrate agarose temperature.
- 4) Once agarose is at 37°C, add 25 µl of semen solution to the vial and mix thoroughly.
- 5) Place a drop of the cell suspension in vial onto the treated face of the slide (marked surface) and cover with a glass coverslip. Avoid making air bubbles by gently pressing. We recommend a drop of 15, 25 or 50 microlitres to coverslips of 18x18 mm, 22x22 mm or 24x60 mm, respectively. Slides need to be always placed in horizontal position.
- 6) Place the slide on the cooled plate into the fridge and leave the sample to solidify for 5 minutes.

Processing the sperm sample

- 7) Remove the coverslip smoothly, and set the slide (in horizontal position) in 10 ml of the lysing solution. This is the volume needed for one slide. Scale the volume depending on the number of slides to be processed. Incubate for 5 min at room temperature (22°C). **Work under fume hood or cover the tray during incubation!**
- 8) Transfer the slide to another plastic tray with distilled water. Wash it for 5 min.
- 9) Dehydrate in sequential 70, 90 and 100% ethanol baths (2 min each) and air dry. Remember to always keep the slide in horizontal position.

10) Once the slides are dried, they can be stored for several months at room temperature in a dry place.

Staining the sample

Proceed to stain the slide just before visualization and analysis. For fluorescence microscopy any DNA-specific binding fluorochrome could be used (DAPI, Propidium Iodide, Ethidium Bromide or Synergy Brand related fluorochromes).

Concentrations of use are identical to those for routine fluorescence microscopy. (i.e. propidium iodide 2,5 µg/ml, and mounted in Vectashield antifading).

Once obtained the concentration, mix the fluorochrome with an antifading solution in proportion 1:1. A final volume of 5 µl is recommended for staining a 25x 25 coverslip.

UNDER THE MICROSCOPE AND CLASSIFYING THE NUCLEOIDS

Any magnification can be used, from 10x to 100x. 20x, 40x or 60x dry objectives are recommended. The nucleoid corresponding to the intensely protein depleted nuclei is made up of two parts: the core, which appears at central position, and the peripheral halo corresponding to the chromatin/DNA dispersion ringlets. The analysis of a minimum number of 500 spermatozooids per sample is recommended. The criteria for classification are as follows:

SPERMATOOZA WITH FRAGMENTED DNA: Spermatozoa with a large and spotty halo of chromatin dispersion.

SPERMATOOZA WITH UNFRAGMENTED DNA: Spermatozoa with a small and compact halo of chromatin dispersion.

SAFETY WARNINGS AND ENVIRONMENTAL PRECAUTIONS

- ☐☐ Biological samples have to be handled as potentially infectious.
- ☐☐ Care should be taken to avoid contact with skin or eyes, and to prevent inhalation. Gloves should be worn to handle the products. Lysing solution contains β-mercaptoethanol. Work under air removal environment and please follow the manufacturer's Material Safety Data Sheet regarding safe handling.
- ☐☐ Do not dispose waste products into the environment.
- ☐☐ Please follow the specific safety regulation of your Faculty or Research Center with respect to chemicals storage and toxic products disposal.


REAGENTS STABILITY AND STORAGE

The components of the kit are stable at room temperature; always keep in a dry and light protected place. Expiration: the reagents supplied are stable for a minimum period of one year. Keep fluorochrome at 4°C.

TROUBLESHOOTING THE PROBLEM -IT'S COMMON CAUSE - OUR SUGGESTED APPROACH

- 1) *Weakly stained halos* - Too short colouring exposure time - Extend time exposure to dye.
- 2) *Failure to differentiate halo and core* - Too long colouring exposure time - Discolour by intense buffer washing and colour again but reducing time exposure to dye.
- 3) *The halos appear displaced from the core in one-way direction* - Slides were not kept in horizontal position during sample processing.

Figure: Arrows show those with fragmented DNA

sperm classification


Removal of sperm nuclear proteins results in nucleoids with a central core and a peripheral halo of dispersed DNA loops. The sperm tails remain preserved. For sperm classification we must take into account that sperm DNA fragmentation is a continuous process which produces a series of different halo sizes. Initially, the study of a minimum of 500 spermatozoa per sample is recommended, adopting the criteria of Fernández et al. (Fertil Steril 84: 833-842, 2005). Avoid scoring cells close to the edge of the microgel. Classification:

- **Spermatozoa without dna fragmentation:**

- **spermatozoa with big halo:** those whose halo width is similar or higher than the minor diameter of the core (Figure 1).
- **spermatozoa with medium-sized halo:** their halo size is between those with large and with very small halo (Figure 2).

- **Spermatozoa with fragmented dna:**

- **spermatozoa with small halo:** the halo width is similar or smaller than 1/3 of the minor diameter of the core (Figure 3).
- **spermatozoa without halo:** (Figure 4).
- **spermatozoa without halo and degraded:** those that show no halo and present a core irregularly or weakly stained (Figure 5).
- **“others”:** cell nuclei which do not correspond to spermatozoa. One of the morphological characteristics which distinguish them is the absence of tail.


Dyn-Halosperm® distinguishes between fragmented and non-fragmented DNA in sperm cells and allows the calculation of the Sperm DNA Fragmentation Index (SDF), i.e. the percentage of spermatozoa with fragmented DNA. In addition, using the Dyn-Halosperm® kit it is possible to study sperm DNA fragmentation over four different time points.

Sperm DNA fragmentation is related to fertilization rate¹, embryo quality¹, embryo development and male pathologies (such as varicocele² or infections by *Chlamydia trachomatis*³).

Each individual has a unique dynamic profile of sperm DNA fragmentation⁴.

Sperm DNA Fragmentation index (SDF):

Information provided by Dyn-Halo sperm*




The Dyn-Halo sperm* protocol provides the following information about DNA fragmentation for a given semen sample:


1. Basal fragmentation level.
2. Rate of change.
3. Mean.

Clinical applications based on basal fragmentation and rate of change

Knowing the evolution of sperm DNA fragmentation over time allows clinicians to:

1. Select the optimum moment to carry out an ART cycle.
2. Select the type of assisted reproduction technique: if the rate of DNA fragmentation is high then it is advisable to use techniques that assure a quick fertilization of the oocyte.
3. Assess the quality of semen samples or donors for suitability.
4. Provide answers to cases of unexplained infertility, ART failure or repeated abortions.

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